


Log In

Go to <http://www.corestudycast.com/>.


Click Client Login in the top right.

The first time you log in, you will be required to change your password and acknowledge the HIPAA policy.



Studies Page

Each time you log in, the Studies page is displayed. This page lists all studies you are authorized to view. As a reading physician, you can view studies that are Preliminary or Final and that are assigned to your reading group.



The screenshot shows the Studycast interface. At the top, there are navigation tabs: 'Studies' (selected), 'Patients', 'Status', 'Admin', 'Help', and 'Support'. On the right, it says 'Logged in as csiadmin-csi | Profile | Logout'. Below the navigation is a search bar and a table of studies. The table has columns for 'Study ID', 'Date of Service', 'Study Type', 'Status', 'Patient Name', 'DOB', and 'Pati'. There are four rows of study data. At the bottom of the table, there are buttons for '+New Study...', 'Save list...', and 'Refresh'. The page number is 'Page 3 of 10' and the view count is 'View 201 - 300 of 1,000'. The footer contains version information and copyright notice.

	Study ID	Date of Service	Study Type	Status	Patient Name	DOB	Pati
<input type="checkbox"/>	60    11298946	2021-04-06 08:58:44	Echo	Preliminary	Zimmerman, Scott	1963-11-13	MHRC-040
<input type="checkbox"/>	56    11298943	2021-04-06 11:49:17	Echo	Preliminary	Fiore, Cynthia	1948-02-09	MHRC-040
<input type="checkbox"/>	45    11298940	2021-04-06 10:50:56	Stress Echo	Final	Carter, Alicia	1958-09-24	MHRC-040
<input type="checkbox"/>	67    11298938	2021-04-06 09:54:41	Echo	Final	Lee, Harrison	1941-12-06	MHRC-040

 To add, remove, or reorder the columns on the Studies page, click the Grid Setup button in the bottom right.

Studies are color coded by status



Preliminary Ready for the reading physician to review.



Final Approved by the reading physician, and the final report has been generated.

The Study ID column displays three clickable icons



Lists the number of images in the study. Click to view the images.




Click to open the worksheet.



Click to view the report (Icon only appears if a report has been created).



Tip If you can't find an exam, talk with the technologist to make sure it's been uploaded, has a status of Preliminary, and has been assigned to your reading group



To filter the list of studies

1. Click **Select Criteria** and select an item from the list.
2. Enter the requested information.
3. Click the **Search** button.

To clear all search criteria and return to the complete list of studies, open the drop-down menu on the **Search** button and select **Reset**.

Complete the Worksheet and Generate the Final Report

Change observations as needed.

To generate findings statements, click the Generate All Findings button.

Click the Generate Conclusions button to generate conclusions based on observations.

To edit findings or conclusions statements, click the text and begin typing.

To enter a new statement, click the dots (...) and begin typing.

To view the report, click Report in the top left or Preview Report in the bottom right.

To approve the study, click either Approve or Approve & Next. If a signature image has been uploaded for you, it will be included in the final report. To upload a signature image, contact your admin.

Access the Studycast system from any Internet-connected device.

For more information, click Help on any page, view the Tutorial Videos, or contact Support.

Phone 866.209.3393 x3

Email support@corestudycast.com



www.corestudycast.com

Copyright © 2022 Core Sound Imaging, Inc.
LBL-72-03-208-B Effective 02/27/22