

Improved patient outcomes, team training and collaboration.

Already at your fingertips.



The Studycast system was designed to streamline your daily clinical workflow, but did you know that it can also help with:

- Peer review
- Training
- Collaboration

Seamlessly integrate peer review

The Studycast system makes it easy to create a fair, consistent process for evaluating clinical performance and improving patient outcomes. With flexible tools at your fingertips, streamlined retrospective peer review can be a non-disruptive part of your day - every day.


- 1 Clone an exam and assign it to a second physician.
- 2 The clone appears in that physician's worklist like any other exam.
- 3 When both interpretations are complete, use the Compare tool for a quick line-by-line comparison.

Want to add a scoring system to rate the discrepancies and generate practice-wide metrics? **Studycast Support can help.**

Improve your team's clinical skills

Maybe you have new techs to train on your protocols and expectations for the preliminary report. Or maybe you have a great case you want to flag for team training and review. Either way, your Studycast account can help.

- Flag model studies that demonstrate best practices by saving a copy to a separate training division.
- Use your flagged studies for monthly or quarterly team training.
- To train a new tech, create a clone of a model study and assign it to the tech. Then, compare the workup with your model line by line using the Compare tool.



The screenshot shows the Studycast web interface. At the top, there is a navigation bar with links for Studies, Patients, Upload Status, Downloads, Admin, Help, and Support. The user is logged in as 'stefanidemo'. Below the navigation bar, there is a search bar and a table of study records. The table has columns for Study ID, Original study ID, Study Type, DOB, Gain, Required Images, Tech Suff, Interp Accurate, Result, and Image grade. There are five rows of data, all for Obstetric studies. The bottom of the interface shows a footer with version information and copyright notices.

Study ID	Original study ID	Study Type	DOB	Gain	Required Images	Tech Suff	Interp Accurate	Result	Image grade
38	9792335	Obstetric	1995-10-10	correct	Yes	No	Yes	tech limited	3
31	9781216	Obstetric	1988-12-28	correct	Yes	No	Yes	tech limited	3
38	9777680	Obstetric	2001-01-16	correct	Yes	No	Yes	tech limited	3
38	9777052	Obstetric	1995-04-13	correct	Yes	Yes	Yes	true positive	4
38	9777005	Obstetric	1987-12-09	correct	Yes	Yes	Yes	true positive	4

Collaborate more effectively

The Studycast system makes it easy to collaborate within your practice and with others outside your facility. In addition to peer review and training, you can

- Flag images and loops of interest for the reading physician
- Leverage remote access to request or provide a second opinion
- Securely share studies with patients or colleagues outside your practice

Want to know more?

Contact Studycast Training at 919.277.0636, option 5



www.corestudycast.com

Studycast is a registered trademark of Core Sound Imaging, Inc.

LBL-72-03-505-A Effective 02/03/2021

© 2021 Core Sound Imaging, Inc. All rights reserved.

5510 Six Forks Road, Suite 200
Raleigh, NC 27609

Office: 919.277.0636

CORESOUND
IMAGING